

ACCURATE
EFFECTIVE
CONVENIENT
FAST

Yikon
GENOMICS · 亿康基因




SARS-CoV-2 Antibody Test (Colloidal Gold)



- An auxiliary and/or supplementary screening kit for detection of SARS-COV-2 Antibodies
- Easy to use, instant result in 15 minutes, Please do not judge the result after 20 minutes.

Antibodies are created upon viral infection. Immunoglobulin M (IgM) is created first, and is therefore an early sign of infection. Immunoglobulin G (IgG) is created later, providing more specific and stronger evidence of infection.



Antibodies can be detected as early as one week after infection, with slight variation among individuals.

1. A negative result indicates no infection, or a very early stage infection
2. A positive result indicates an infection with coronavirus.

Order information

Product name	Packing size	Sample type	Storage condition
SARS-CoV-2Antibody Test (Lateral flow method)	1 test per kit 10 tests per kit 100 tests per kit	20 shares per kit 50 shares per kit	Whole blood, serum and plasma 2°C~ 30°C

Yikon
GENOMICS · 亿康基因

Phone: (800) 674-7781
Email: ussales@yikongenomics.com
brandon@yikongenomics.com
Address: 16192 Coastal Hwy Lewes DE 19958

